

RNA prep (12/20/10)

The protocol for harvesting, preparing, and labeling RNA was taken from Maitreya Dunham's website

Hybrid of the DeRisi protocol (www.microarrays.org) and a standard acid-phenol prep circulating around the Brown/Botstein labs circa 2001

Use RNase free reagents and glass-/plastic-ware throughout! Remember to use glass pipets with chloroform.

Lysis buffer (100 ml)

2 ml	0.5 M EDTA
5 ml	10% SDS
1 ml	1 M Tris pH 7.5
92 ml	RNase-free water

Remove a manageable set of samples from the -80. They should be in 2 ml eppendorf tubes.

Before they thaw, add 750 μ l lysis buffer. Vortex, trying to get all the cells off the membrane.

Add 750 μ l acid phenol. Vortex.

Incubate 1 hour 65C, vortexing every 20 minutes.

Fish out the filter and discard.

Ice 10 min.

While they are incubating, spin the 2 ml heavy phase lock gel tubes (made by Eppendorf, sold lots of places) for 30 sec full speed in a room temperature microcentrifuge. Set aside.

Spin lysate 5 min.

With a pipet, transfer the top aqueous layer to the PLG tube.

Add 750 μ l chloroform. Invert to mix. Do not vortex!

Spin 5 min.

Pour aqueous layer into a new 15 ml falcon tube.

Add 75 μ l (or 1/10 volume if you lost some) 3 M sodium acetate. Mix.

Add 1.5 ml (or 2 volumes) ethanol. Mix.

Incubate -20C >30 min (preferably overnight).

Spin 3000 rpm 10 min.

Wash pellet 2X with 70% ethanol, with 2 min 3000 rpm spins between washes.

Air dry inverted on the bench 30 min.

Dissolve pellet in 100 μ l water.

Measure concentration of crude RNA with nanodrop.

Clean up the RNA prep with QIAGEN RNeasy kit. Elute in 30 μ l water.

Measure concentration of clean RNA with nanodrop.

Optional: Check the quality of the RNA on the Bioanalyzer (the Bioanalyzer requires 2 μ l of 50 – 500 ng/ μ l).