

PCR (12/23/08)

For a 50 ul Taq PCR:

40 ul H₂O
5 ul 10 x Taq buffer
1.5 ul dNTPs (10 uM each dNTP)
1 ul Primer 1 (usually 20 – 40 uM)
1 ul Primer 2 (usually 20 – 40 uM)
1 ul Template
0.5 ul Taq Polymerase
50 ul TOTAL

Thermal Cycler Protocol:

1 96° C 2 min *Hot Start*
2 96° C 30 s *Denaturing*
50° C 30 s *Annealing*
72° C 2 min *Extension*
cycle 30 x
3 72° C 5 min *Final extension*

Variations on the thermal cycler protocol: For any given PCR, I try the standard protocol first. Typically, I will troubleshoot a failed PCR by modifying the annealing temperature or the extension time. The rule of thumb for the annealing temperature is the $T_m - 5^\circ \text{C}$. The rule of thumb for the extension time is 1 min for every kilobase to be amplified.

10 x Taq Buffer

500 mM KCl
10 mM Tris-Cl, pH 8.4
1 mg/ml Gelatin
25 mM MgCl₂
0.5% Tween 20