

Yeast genomic preparation for sequencing (08/30/17)

This protocol starts with cell pellets from 5 ml of culture.

Resuspend cells in 300 μ l of Lysis buffer.

Add ~0.3 g Glass beads

Add 300 μ l of phenol:chloroform:isoamyl alcohol

Vortex for 10 min.

Spin 5 min at top speed.

Remove aqueous layer (~250 μ l) and transfer to a new tube.

Add 5 μ l RNase (100 mg/ml). Incubate 37°C for at least 1 hour.

Add 300 μ L of phenol:chloroform:isoamyl alcohol. Vortex for 15 seconds.

Transfer contents to a phase-lock gel (light) tube.

Spin 5 min at top speed.

Dump aqueous layer into a tube with 500 μ l of ethanol. Precipitate DNA in -20° C freezer for at least 30 min (overnight).

Pellet DNA by spinning for 10 minutes at top speed. Decant or aspirate ethanol.

Add 1 ml of 70% ethanol, vortex, and spin for 2 min at top speed. Decant or aspirate ethanol.

Dry DNA pellet in speed-vac for ~10 min. Resuspend pellet in 100 μ l of ddH₂O.

Check and record DNA concentration using Qubit fluorometer

Lysis Buffer

2% Triton X-100

1% SDS

100 mM NaCl

10 mM Tris-Cl, pH 8.0

1 mM EDTA

To make 50 ml of Lysis Buffer

40.9 ml H₂O

5 ml of 20% Triton X-100

2.5 ml of 20% SDS

1 ml of 5 M NaCl

500 μ l of 1 M Tris-Cl, pH 8.0

100 μ l of 0.5 M EDTA