

RNA labeling (12/5/07)

The protocol for harvesting, preparing, and labeling RNA was taken from Maitreya Dunham's website

The amplification/labeling is done per Agilent instructions with half volume reactions and a quarter amount of recommended dye.

Make a 100 ng/μl stock of total RNA.

Aliquot 3.25 μl (325 ng) total RNA into a PCR tube.

Add 2.5 μl T7 Promoter Primer.

65C 10 min

ice 5 min

Warm 5X first strand buffer at 100C, with occasional vortexing, until it completely dissolves (~1 min).

Prepare cDNA master mix, in this order at RT:

1X	
2 μl	5X first strand buffer
1 μl	0.1 M DTT
0.5 μl	10 mM dNTP mix
0.5 μl	MMLV RT
0.25 μl	RNaseOUT

Add 4.25 μl to each reaction.

40C 2 hours

65C 15 min

ice 5 min

Add 0.8 μl Cy3 CTP to each reference reaction.

Add 0.6 μl Cy5 CTP to each experimental reaction.

Warm 50% PEG at 40C until it's resuspended and easy to pipet.

Prepare transcription master mix, in this order at RT:

1X	
8.25 μl	water
10 μl	4X transcription buffer
3 μl	0.1 M DTT
4 μl	NTP mix
3.2 μl	50% PEG
0.25 μl	RNaseOUT
0.3 μl	inorganic pyrophosphatase
0.4 μl	T7 RNA polymerase

Mix by pipetting.

Add 29.4 μl to each reaction. Mix by pipetting.

40C 2 hours in the dark

Purify with an RNeasy column.

Nanodrop 1 ul to check yield and dye incorporation.

Store RNA at -20C.